dPOFLEX Industrial Peristaltic Pump GP01, BP01



Pump drive model	dPOFLEX GP01	dPOFLEX BP01			
Housing material	Cast aluminum	304sst			
Speed range	0.1rpm-265rpm, CW/CCW				
Speed resolution	0.1rpm				
Display	4.3 inch LCD, Chinese or English				
Clock function	Year/month/day/hour/minute can be set, hours and minutes are displayed on the running screen				
Work mode	Flow rate mode, dispensing mode, calibration mode				
Flow rate mode	Pump runs continuously at the set flow rate, and displays flow rate and pump speed on the screen				
Flow rate setting range	1.1ml/min-17.0L/min (depending on pump head and tubing)				
Dispensing mode	Pump repeats the dispensing at set volume, cycles and interval time, and displays single volume, dispensed volume, dispensed cycles/target cycles on screen				
Volume of single dispense	1.0mL-274.9L (depending on pump head and tubing)				
Dispensing cycles	1-9999 (0 for unlimited)				
Interval time	5s-9999s				
Calibration function	Improve the flow rate and dispensing volume accuracy				
Control mode	Keypad control, external signal control and communication control				
Communication control	Rs485 or TCP/IP(optional) interface, Modbus RTU, Industrial Ethernet(optional), Profibus DP (optional)				
External control	Speed control: 0-10V, 4-20mA Start/stop, direction control: 5V-24V logic level signal or relay signal (dry contact)				
Status output	Pump Speed: 0-10V, 4-20mA, frequency signal (optional) Running/stop, direction, tubing leak alarm: relay signal (dry contact)				
Open head sensor	Pump will stop when pump head is open for enhanced user safety				
Status when power up	Stop				
Keypad lock function	Enable the keypad lock function to lock the keypad automatically after 20s standby, (start/stop key will not be locked)				
Tubing leak detection	Detect the tubing leak and give alarm (optional)				
Flow rate measurement	Measure the flow rate by connecting flow meter (optional)				
User access control	Three levels of user accesses (administrator, technician, operator)				
Electronic signature and audit trail	Meet the requirement of FDA 21CFR Part 11				
Log record	Record 90 days operation history, and can be exported directly				
Support IQ/OQ	Meet GMP requirement				
Certificate	CE, UKCA EN 61010-1:2010, EN IEC 61010-2-201; EMC: EN IEC 61000-6-2, EN IEC 61000-6-4				
Deminsion(LxWxH)	260x407x298mm	260x407x290mm			
Power supply	AC100-240V 50/60Hz				

≤17.6kg

Max. Flow Rate:



🕨 17L/min

Available in stainless steel housing and cast aluminum housing for variety of applications

IP66 for wet and dusty environment

Supporting fieldbus communication control, digital signal control, analog signal control, and pump status output

Three-level user authority management, with electronic signature and audit trail functions, meet 21CFR Part 11 and GMP lab requirements

Integrating multiple sensors: fluid leakage sensor, open-head sensor and flow rate sensor, to realize predictive maintenance

BP01+LPH01





BP01+GPH01





GP01+LPH01









Model (Product Code)	Pump Head	Roller Number	Tubing	Reference Max. Flow Rate (L/min)
	GPH01	2	26#,73#,82#,184#	15
dPOFLEX GP01 (05.02.80A)	GPH02	2	186#,188#	17
dPOFLEX BP01 (05.02.81A)	GPH03	4	186#,188#	14
	LPH01	6	185#,186#,187#	13.5

Temperature: 5°C-40°C, relative humidity≤100% IP 66

Operating condition

IP rating Weight (drive only)



≤16.6kg





E - mail:info@longerpump.com

13.5 (double tube

assembly)

Fax: 0312 - 3168553 Internet

187#(16x4)