

GUARANTEED TO IMPRESS

Accumax Electronic Pipettes



User-friendly
Graphical
Interface



3 Point Precise
Auto-calibration



5 Different
Pipetting Modes

- Forward pipetting
- Reverse pipetting
- Multi-Dispensing
- Diluting • Mixing




Light Weight
Construction



Calibrated As Per ISO
8655 In ISO 17045
Accredited Lab



Ultra-efficient
Battery



Individual Auto-sensing
Magnetic Piston For Each
Channel Ensures Precise Pipetting



Intelligent
Sleep Mode



* Also Available In Multichannel

For more details visit us at www.accumaximum.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.